

Standard Light Source for Operation and Sensitivity Check

Lumitester Calibration Kit

Smart/PD-30/PD-20

Product Code : 60484

※This device is factory-calibrated, with a minimum period of validity as stated on the label.
 ※Following the stated period of validity, the device can be re-calibrated and re-labeled with new period of validity for a fee.



Positive Rod

Negative Rod

Lumitester Smart

Specification

Light Source - Positive Rod	
Light emitting device	Green LED
RLU value	Displayed on the label of light source, (positive rod)
Charging time	30 seconds or more
Dimensions	Diameter: ϕ 12mm Length: 135mm
Weight	About 20g


Blank Tube – Negative Rod	
RLU value	Typ. 0RLU
Tolerance of RLU value	Maximum 10RLU
Dimensions	Diameter: ϕ 12mm Length: 150mm
Weight	About 10g

※ This equipment is Lumitester Smart/PD-30/PD-20 exclusive use.
 ※ Please refer to the instruction manual for full procedures.
 (Instruction manual is available in the website.)

Accessories



Usage (positive rod)



Kikkoman Biochemifa Company

2-1-1 Nishi-shinbashi Minato-Ku, Tokyo, 105-0003

E-mail : biochemifa@mail.kikkoman.co.jp
 URL : <https://biochemifa.kikkoman.com/e/>

※The contents described in this leaflet are subject to change without notice.
 ※“Lumitester” is the registered trademark of Kikkoman Corporation in Japan and other countries.

© 2020 Kikkoman Corp. (2502C230101)